

CELSIUS® Ablation Catheters Featuring Thermocouple (TC)

CELSIUS® THERMOCOOL® Catheters




2-5-2 mm

10-Pin Connector

Associated Cable:
Catalog No: 39E43R (Stockert EP Shuttle Generator)



Associated Tubing:
Catalog No: CFT-001 (CoolFlow™ Tubing)



European Ordering No.	F	Curve	Tip Electrode	Spacing	Insertion Length	Color Code	Stockert EP Shuttle Interface Cable
35Q13R	7F	B	3.5 mm	2-5-2 mm	115 cm	red	39E43R
35Q33R	7F	D	3.5 mm	2-5-2 mm	115 cm	blue	39E43R
35Q53R	7F	F	3.5 mm	2-5-2 mm	115 cm	orange	39E43R
35Q73R*	7F	J	3.5 mm	2-5-2 mm	115 cm	black	39E43R

* For availability, please contact your local Biosense Webster representative.

EZ STEER® THERMOCOOL® Catheters Featuring Thermocouple (TC)

EZ STEER® THERMOCOOL® Catheters




10-Pin Connector

Associated Cable:
Catalog No: 39E43R (Stockert EP Shuttle Generator)



Associated Tubing:
Catalog No: CFT-001 (CoolFlow™ Tubing)



European Ordering No.	F	Curve	Tip Electrode	Spacing	Insertion Length	Color Code	Stockert EP Shuttle Interface Cable
36Q33M	7.5F	D-D	3.5 mm	2-5-2 mm	115 cm	blue/blue	39E43R
36Q35M	7.5F	D-F	3.5 mm	2-5-2 mm	115 cm	blue/orange	39E43R
36Q55M	7.5F	F-F	3.5 mm	2-5-2 mm	115 cm	orange/orange	39E43R
36Q5JM	7.5F	F-J	3.5 mm	2-5-2 mm	115 cm	orange/black	39E43R
36QJIM	7.5F	J-J	3.5 mm	2-5-2 mm	115 cm	black/black	39E43R

NAVISTAR® Ablation Catheters featuring Thermocouple (TC)

NAVISTAR® THERMOCOOL® Catheters



2-5-2 mm

25-Pin Connector



Associated Cable: Catheter to CARTO® 3 System PIU
Catalog No: CR3425CT



To Catheter
Black 25-pin


To PIU
Red 34-pin



CARTO® 3 System PIU to Stockert EP Shuttle Generator
Catalog No: 39E43R



Associated Tubing:
Catalog No: CFT-001 (CoolFlow™ Tubing)



European Ordering No.	F	Curve	Tip Electrode	Spacing	Insertion Length	Color Code	Interface Cable
34H17M	7.5F ²	B	3.5 mm	2-5-2 mm	115 cm	red	CR3425CT
34H27M	7.5F ²	C	3.5 mm	2-5-2 mm	115 cm	green	CR3425CT
34H37M	7.5F ²	D	3.5 mm	2-5-2 mm	115 cm	blue	CR3425CT
34H57M	7.5F ²	F	3.5 mm	2-5-2 mm	115 cm	orange	CR3425CT
34HJ7M*	7.5F ²	J	3.5 mm	2-5-2 mm	115 cm	black	CR3425CT

² 8 french ring and tip electrodes.

* For availability, please contact your local Biosense Webster representative.

EZ STEER® THERMOCOOL® Nav Catheters Featuring Thermocouple (TC)

EZ STEER® THERMOCOOL® NAV Catheters



2-5-2 mm

25-Pin Connector

Associated Cable: Catheter to CARTO® 3 System PIU
Catalog No: CR3425CT



Note: With 2 additional shaft electrodes



To Catheter
Black 25-pin


To PIU
Red 34-pin



CARTO® 3 System PIU to Stockert EP Shuttle Generator
Catalog No: 39E43R



Associated Tubing:
Catalog No: CFT-001 (CoolFlow™ Tubing)



European Ordering No.	F	Curve	Tip Electrode	Spacing	Insertion Length	Color Code	Interface Cable
36H33M	7.5F	D-D	3.5 mm	2-5-2 mm	115 cm	blue/blue	CR3425CT
36H35M	7.5F	D-F	3.5 mm	2-5-2 mm	115 cm	blue/orange	CR3425CT
36H55M	7.5F	F-F	3.5 mm	2-5-2 mm	115 cm	orange/orange	CR3425CT
36H5JM	7.5F	F-J	3.5 mm	2-5-2 mm	115 cm	orange/black	CR3425CT
36HJJM	7.5F	J-J	3.5 mm	2-5-2 mm	115 cm	black/black	CR3425CT

NAVISTAR® Ablation Catheters Featuring Thermocouple (TC)

NAVISTAR® THERMOCOOL® Catheters



2-5-2 mm

➔ 25-Pin Connector


Associated Cable: Catheter to CARTO® XP System PIU
Catalog No: 39E58M



CARTO® XP System PIU to Stockert EP Shuttle Generator
Catalog No: 39E43R



Associated Tubing:
Catalog No: CFT-001 (CoolFlow™ Tubing)



European Ordering No.	F	Curve	Tip Electrode	Spacing	Insertion Length	Color Code	Interface Cable
34H17M	7.5F ²	B	3.5 mm	2-5-2 mm	115 cm	red	39E58M
34H27M	7.5F ²	C	3.5 mm	2-5-2 mm	115 cm	green	39E58M
34H37M	7.5F ²	D	3.5 mm	2-5-2 mm	115 cm	blue	39E58M
34H57M	7.5F ²	F	3.5 mm	2-5-2 mm	115 cm	orange	39E58M
34HJ7M*	7.5F ²	J	3.5 mm	2-5-2 mm	115 cm	black	39E58M

² 8 french ring and tip electrodes.

* For availability, please contact your local Biosense Webster representative.

EZ STEER® THERMOCOOL® Nav Catheters Featuring Thermocouple (TC)

EZ STEER® THERMOCOOL® Nav Catheters



2-5-2 mm

25-Pin Connector


Associated Cable: Catheter to CARTO® XP System PIU
Catalog No: 39E58M



CARTO® XP System PIU to Stockert EP Shuttle Generator
Catalog No: 39E43R



Associated Tubing:
Catalog No: CFT-001 (CoolFlow™ Tubing)



European Ordering No.	F	Curve	Tip Electrode	Spacing	Insertion Length	Color Code	Interface Cable
36H33M	7.5F	D-D	3.5 mm	2-5-2 mm	115 cm	blue/blue	39E58M
36H35M	7.5F	D-F	3.5 mm	2-5-2 mm	115 cm	blue/orange	39E58M
36H55M	7.5F	F-F	3.5 mm	2-5-2 mm	115 cm	orange/orange	39E58M
36H5JM	7.5F	F-J	3.5 mm	2-5-2 mm	115 cm	orange/black	39E58M
36HJJM	7.5F	J-J	3.5 mm	2-5-2 mm	115 cm	black/black	39E58M

THERMOCOOL® SF Bi-Directional NAV Catheters Featuring Thermocouple (TC)

THERMOCOOL® SF Bi-Directional NAV Catheter*

➔ 25-Pin Connector


Associated Cable: Catheter to CARTO® XP System PIU
Catalog No: 39E58M



CARTO® XP System PIU to Stockert EP Shuttle Generator
Catalog No: 39E43R



Associated Tubing:
Catalog No: CFT-001 (CoolFlow™ Tubing)



European Ordering No.	Curve Type	French size (F)	Number of Electrodes	Spacing	Tip Electrode	Insertion Length	Interface Cable
BNI35BBH	BB	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35BDH	BD	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35BFH	BF	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35DDH	DD	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35DFH	DF	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35DJH	DJ	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35FFH	FF	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35FJH	FJ	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M
BNI35JJH	JJ	8	4 + 2	2-5-2 mm	3.5 mm	115 cm	39E58M

* THERMOCOOL® SF requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.