



CRISTACATH®, HALO® XP, & ISMUS™ CATH Catheters Twenty Electrode with Redel 10-Pin Connector



(2) Redel 10-Pin Connector

Associated Cable:
Catalog No: (2) C6MRMST10SA
European Order No: (2) 39F-54R or 39F-41R



New Catalog No.	Made to Order	Old Catalog No.	Tip Electrode	Size	Catheter Curve	Spacing	Insertion Length	European Order No.
D7A20131RT		D7-A20-131-RT	3 mm	7F	D-Type/CRISTACATH®	1-3-1 mm*	115 cm	36Y-39I
D7R20P14RT		D7R-20-P14-RT	None	7F	13 cm/ISMUS™ CATH	2-12-2 mm	115 cm	36Y-32I
D7T20282RT		D7-T20-282-RT	None	7F	Tricuspid/HALO® XP 20 mm Tail	(20) - 2-8-2 mm	110 cm	36J-13I
D7T20P15RT		D7-T20-P15-RT	None	7F	Tricuspid/HALO® XP	2-13-2 mm	110 cm	36J-16I
116008RT		1160-08-RT	None	7F	Tricuspid/HALO® XP 20 mm Tail	1-18-2-8-2 mm	110 cm	36J-14I