

Perclose ProGlide

Suture-Mediated Closure System

Simple

Secure¹

Control

The security¹
of suture



Simple

- Pre-tied knot
- Numbered deployment sequence
- Single operator procedure with QuickCut mechanism
- 5-8F

Secure¹

- Secure¹ close with a polypropylene monofilament suture
- Minimized inflammatory response
- The closure can be challenged and confirmed on the table
- No reaccess restrictions if previous arteriotomy repairs were achieved with Perclose ProGlide

Control

- Ability to maintain wire access
- Reduced time to hemostasis, ambulation and discharge



**Endovascular
Commitment**

Abbott
A Promise for Life



The security¹ of suture



Pre-tied knot



Single operator procedure with QuickCut mechanism



Numbered deployment sequence



Secure¹ close with a polypropylene monofilament suture




Guide wire exit port allows wire access to be maintained



Foot deployment



Needles deployment



Connection of suture with pre-tied knot



Suture brings the tissues together to promote primary healing – no need for a clot to form to close the arteriotomy

Ordering Information – please fill in your required quantity and return to your Abbott representative

Ref. No.	Quantity	Description	Quantity
12673-05	<input type="text"/>	Includes: • 1 Perclose ProGlide 6F Suture-Mediated Closure Device • 1 Suture Trimmer	10 systems per box

Representative name	<input type="text"/>	Contact number	<input type="text"/>	Physician's signature	<input type="text"/>
---------------------	----------------------	----------------	----------------------	-----------------------	----------------------

Reference: 1. Source: Data on file at Abbott Vascular

All products herein are trademarks of the Abbott Group of Companies. Perclose ProGlide is subject to prior training requirement as per the Instructions for Use.

Contact your local representative for more information

Photo(s) on file at Abbott Vascular. All drawings are artist's representations only and should not be considered as an engineering drawing or photograph. Products intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for instructions for use, warnings and potential complications associated with use of this device. Information contained herein is for distribution for Europe, Middle East and Africa ONLY. Please check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force.

For more information, visit our website at abbottvascular.com

© 2012 Abbott. All rights reserved. 2-EH-1-0459-01 04/2012 Rev. C

