



Armada 14 XT

PTA Catheter

Redefining possibilities
in BTK therapy



Ability to perform in challenging situations

Dual Hydrophilic Coating

Two layers of outer coating for outstanding durability and to maintain lubricity during long, complex cases

Low Crossing Profile

Excellent for pre-dilatation and designed for hard-to-cross lesions



Transitionless Tip

For confident crossing in resistant lesions




Armada 14 XT

PTA Catheter


Ability to perform in challenging situations

89% better overall performance in extremely and moderately challenging procedures


Overall Performance



Procedural Challenge



Source: Physician Preference Testing RPT2091726 Rev A. Data on file at Abbott Vascular. Overall Performance as compared to Amphirion Deep, Apex, Sterling ES, Coyote, Falcon, Nanocross, Sleek, Journey, Pacific XT, Advance 14 LP and Senri.



Ordering Information

	Balloon Diameter (mm)	Balloon Length (mm)	Stock Number (Catheter Length 145 cm)
Armada 14 XT OTW	1.2	12	F2012-012
	1.2	20	F2012-020
	1.5	12	F2015-012
	1.5	20	F2015-020
	2.0	12	F2020-012
	2.0	20	F2020-020

Technical Information	
Configuration	OTW
Guide Wire Compatibility	0.014"
Sheath Compatibility	4F
Nominal Pressure	8 atm
Rated Burst Pressure	14 atm

Armada 14 XT is a trademark of the Abbott Group of Companies.

Apex, Sterling ES and Coyote are trademarks of Boston Scientific Corporation or its affiliates.

Sleek and Journey are trademarks of Cordis Corporation, a Johnson & Johnson Company.

Advance 14 LP is trademark of Cook, Inc.

Nanocross is a trademark of Covidien, Inc.

Amphirion Deep, Pacific XT and Falcon are trademarks of Medtronic, Inc.

Senri is trademark of Terumo.

Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

Photo(s) on file at Abbott Vascular. All drawings are artist's representations only and should not be considered as an engineering drawing or photograph. Products intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for instructions for use, warnings and potential complications associated with use of this device. Information contained herein is for distribution for Europe, Middle East and Africa ONLY. Please check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force.

For more information, visit our website at abbottvascular.com
© 2013 Abbott. All rights reserved. 2-EH-1-3100-01 01/2013 Rev. A