

DATASHEET

Because of the innovative quality design of the Ambu® King Mask, it is the market leading anaesthesia mask design available today. The Ambu® King Mask has a thin and tacky cushion which improves the ability to create a ventilation seal on a variety of patients. The design with rounded edges secures a smooth grip and user comfort.

The Ambu® King Mask is made of phthalate-free PVC material for patient and user safety. Furthermore, the Ambu® King Mask can be used in a MR environment, as it is MR Safe.

The Ambu® King Mask is also available in a scented version with Fresh Scent, Bubble Gum, Strawberry and Cherry scents. The Sweet Dreams™ masks release soothing aroma from natural oils in mask cushion. This can have a calming effect during induction, especially for children¹.

Five Scents

- Unscented
- Fresh Scent (mint)
- Bubble Gum
- Strawberry
- Cherry

Seven Sizes

- Size 1 Neonate
- Size 2 Infant
- Size 3 Toddler
- Size 4 Child / Small Adult
- Size 5 Medium Adult
- Size 6 Large Adult
- Size 7RD XL Adult

Ambu® King Mask

Disposable anaesthesia face mask



KEY BENEFITS

Market Leading Design

Cushion is pliable and tacky to the touch, for an excellent face seal

Ribbed on top of the cushion, for a no-slip grip

Anatomically shaped cuff enabling a tight seal with minimum applied pressure

Flexible dome

Clear dome for easy observation of the patient's condition

Latex free








MR Safe

Phthalate-free material






Available with or without hook ring*

¹ Carson KD et al. Quantitative assessment of patient preference for anaesthesia face mask. Irish Journal of Medical Science 1995;164 suppl. 14
* Not available for sizes 1 and 2.

ORDERING SPECIFICATIONS WITHOUT HOOK RING

Mask	Size 1 Neonate	Size 2 Infant	Size 3 Toddler	Size 4 Child / Small Adult	Size 5 Medium Adult	Size 6 Large Adult	Size 7RD XL Adult
Packaging Color Code							
Original (unscented)	1012	1022	1032 H	1045 H	1055 H	1065 H	1075 H
Fresh Scent (mint)	2012	2022	2032 H	2045 H	2055 H	2065 H	2075 H
Bubble Gum		K1122	1132 H	1145 H			
Strawberry		1222	1232 H	1245 H			
Cherry		1322	1332 H	1345 H			
Scented masks variety pack		1423					

WITH HOOK RING

Mask	Size 3 Toddler	Size 4 Child / Small Adult	Size 5 Medium Adult	Size 6 Large Adult	Size 7RD XL Adult
Packaging Color Code					
Original (unscented)	1032	1045	1055	1065	1075
Fresh Scent (mint)	2032	2045	2055	2065	2075
Bubble Gum	1132	1145			
Strawberry	1232	1245			
Cherry	1332	1345			
Scented masks variety pack	1433	1443			

Part #'s ending in 2 = 20 per case. Part #'s ending in 5 = 50 per case

Part #'s ending in 3 = 30 per case (Scented mask variety pack 10 of each of Bubble Gum, Strawberry, & Cherry)

MATERIALS

Product	Material
Cushion	PVC Polyvinyl Chloride (Phthalate-free)
Crown	Polycarbonate
Hook ring	Polycarbonate
Valve Body	Polycarbonate
Valve Stem	Silicone Rubber
Valve Plug	PP Polypropylene

Packaging

Pouch	LDPE
Box	Cardboard

Environment

Phthalate-free	Yes
Not made with natural rubber latex	Yes

PRODUCT WEIGHT


Size	Weight
1	15 g
2	18.4 g
3	32.4 g
4	38.3 g
5	50.6 g
6	57.9 g
7RD	51.3 g



Unscented size 5



Fresh Scent



Strawberry



Cherry



Bubble Gum

Ambu

 Ambu A/S

Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
ambu.com

LBL-007129 - V07 - 2021/10 - Ambu A/S. Technical data may be modified without further notice. TCC 11173

 2797 US: Rx only